

# BMA101 (DIBMA) User Guide



Despite the wide use of conventional detergents, removal of a membrane protein from the native lipid bilayer and into a detergent micelle (Figure 1) can possibly interfere with the structure, function, and stability of the protein<sup>(1,2)</sup>. To avoid some of these issues with detergents, many systems have been developed to mimic the native lipid membrane. These include bicelles (Figure 2)<sup>(3)</sup>, MSP-nanodiscs (Figure 3)<sup>(4)</sup>, and detergent-free polymer systems such as styrene-maleic acid co-polymers (SMAs)<sup>(5)</sup> and diisobutylene-maleic acid (DIBMA)<sup>(6,7)</sup> (Figure 4). DIBMA from Anatrace is a highly purified version of diisobutylene and maleic acid co-polymer in a 1:1 ratio.




Figure 1




Figure 2




Figure 3




Figure 4

This user guide contains a general protocol and suggestions on how to use the DIBMA polymer to solubilize, purify, and characterize membrane proteins.

## Product Information:

- ★ SKU: BMA101
- ★ Name: Diisobutylene Maleic Acid co-polymer, sodium salt / DIBMA / Purified Sokalan® CP9<sup>(8)</sup>
- ★ MW: ~12,000 g/mol
- ★ dn/dc: 1.346 M<sup>-1</sup>
- ★ Solubility: > 10% in H<sub>2</sub>O
- ★ Absorbance at 280 nm (1% solution): < 0.3
- ★ Mg<sup>2+</sup> Tolerance: 25 mM
- ★ Ca<sup>2+</sup> Tolerance: 20 mM

## DIBMA Protocol:

- ★ To solubilize membrane proteins from the native lipid bilayer, DIBMA can be added directly to a preparation of isolated membranes that have been resuspended in buffer. The DIBMA powder can be added to the membranes to a final concentration of 1% – 3%.
- ★ Alternatively, a stock solution of DIBMA in buffer can be prepared and added to resuspended membranes to a final concentration of 1% - 3%.
- ★ Solubilization can be performed at room temperature or at 4°C with gentle rocking. When compared to solubilization with other polymers, such as SMA, typically a slightly longer solubilization time is needed for optimal extraction efficiency.

# BMA101 (DIBMA) User Guide

- ★ The pH of the solubilization buffer has been shown to effect solubilization efficiency<sup>(6)</sup>, thus small scale solubilization experiments should be performed to determine the optimal pH.
- ★ Upon solubilization, no additional DIBMA polymer needs to be present in buffers used for protein purification, and biophysical analysis.
- ★ Solubilized membrane proteins in DIBMA are compatible with all commonly used protein purification techniques, including IMAC, affinity, ion exchange, and SEC.
- ★ Membrane proteins solubilized and purified in DIBMA are compatible with most biophysical characterization assays, including fSEC, SEC-MALS, UV/Vis spectroscopy, western blotting, and SDS-PAGE (see below).
- ★ Compared to SMA, the absorbance of DIBMA at 280 nm is significantly lower due to the absence of aromatic groups. Additionally, DIBMA is tolerant to higher concentrations of divalent cations.
- ★ The addition of DIBMA to solubilize membrane proteins from whole cells has been performed with limited success.

## SDS-PAGE Protocol:

The presence of polymers (both SMA and DIBMA), may result in the smearing of bands on SDS-PAGE gels<sup>(6,9)</sup>. Proteins can be separated from polymers through a simple precipitation using methanol, chloroform, and water. This protocol was originally described in<sup>(10)</sup>, and slightly modified in<sup>(6)</sup> for use with the DIBMA polymer.

- ★ Measure volume of polymer-containing protein sample
- ★ Add 4x volume of cold methanol to sample and mix thoroughly via vortexing
- ★ Add 1x volume of cold chloroform to sample and mix thoroughly via vortexing
- ★ Add 3x volume of cold water, mix, and centrifuge for 2 min at 4°C and 14,000 g.
- ★ Remove top aqueous layer, add 4x volume of methanol, and mix
- ★ Pellet proteins by centrifuging at 5,000 g for 1 min and 20,000 g for 5 min at 4°C
- ★ Remove supernatant and dry pellet
- ★ Add SDS loading buffer, run SDS-PAGE as normal.

## References:

1. Privé G. G. (2007) *Methods* 41(4), 388-397.
2. Seddon, A. M. et al. (2004) *Biochim Biophys Acta*. 1666(1-2), 105-117.
3. Sanders, C. R., et al. (1994) *Progress in NMR Spectroscopy* 26(5), 421-444.
4. Bayburt, T. H. and Sligar, S., G. (2003) *Protein Sci.* 12(11), 2476-2481.
5. Knowles, T. J. et al. (2009) *J. Am Chem Soc.* 131(22), 7484-7485.
6. Oluwole, A. O. et al. (2017) *Angew Chem Int Ed Engl.* 56(7), 1919-1924.
7. Oluwole, A. O. et al. (2017) *Langmuir* 33(50), 14378-14388.
8. Sokolan is a registered trademark of BASF Corporation
9. Lee, S. C. et al. (2016) *Nature Protocols* 11(7) 1149-1162.
10. Wessel D. and Flügge U. I. (1984) *Anal. Biochem.* 138 141-143.